## SARS-CoV-2 MonlabTest®

**MO-804032 20 TESTS**

*One Step to detect SARS-CoV-2 Antigen*

IVD

A rapid one step test for the qualitative detection of nucleoprotein antigen of SARS-CoV-2 from human nasopharyngeal swab specimens.

For professional *in vitro* diagnostic use only.

### INTENDED USE

The SARS-CoV-2 MonlabTest® is a rapid chromatographic immunoassay for the qualitative detection of Coronavirus (SARS-CoV-2) antigens in human nasopharyngeal swab specimens to aid in the diagnosis of Coronavirus (SARS-CoV-2) infection.

### SYNTHESIS

Recently, SARS-CoV-2 (Severe Acute Respiratory Syndrome Corona Virus 2) emerged as a new specie able to infect humans and currently causing the COVID-19 (Coronavirus disease-19) pandemic, with an overwhelming hospitalization of infected patients.

Clinically, patients with SAR-CoV-2 infection tend to suffer symptoms such as olfactory and gustatory dysfunction, fever, dry cough, anosmia, fatigue, dyspnoea, headache, diarrhoea and sore throat, followed by vascular and systemic complications, such as leukocyte infiltration of the lungs. COVID-19 commonly results in pneumonia, which can evolve to into acute respiratory distress syndrome, leading to respiratory or multiorgan failure.

Coronaviruses (CoVs) belongs to the order of Nidovirales, identified by its envelope characteristics and positive-sense RNA as genetic material. CoVs consist of four structural proteins: spike (S), membrane (M), envelope (E) and nucleocapsid (N) proteins.

The M and E proteins are essential for virus assembly, while the S protein, on the surface of the viral particles, is crucial for affinity and attachment to host cells. The S protein, responsible for viral entry, determines the host tropism and virus transmission. The N protein, the main structural protein of SARS-CoV-2 is responsible for the transcription and replication of viral RNA, the packaging of the encapsulated genome into virions and interactions with the cell cycle of host cells. In addition, the N protein, which has a substantial immunogenic ability, is abundantly expressed during viral infection.

It is believed that the high transmissibility of COVID-19 is related to its high viral loads in the upper respiratory tract and the fact that many individuals remain asymptomatic, shedding and transmitting the virus. The COVID-19 virus transmission occurs via droplets of saliva or discharge from the nose of infected individuals, therefore following proper hygiene practices when coughing and sneezing are keys to deplete transmission.

### PRINCIPLE

The SARS-CoV-2 MonlabTest® is a qualitative lateral flow immunoassay for the detection of Coronavirus (SARS-CoV-2) antigen in human nasopharyngeal samples. The membrane is pre-coated with monoclonal antibodies against SARS-CoV-2 antigens on the test line region. During testing, the sample reacts with the particle coated with anti-SARS-CoV-2 antibodies which was pre-dried on the test strip. The mixture moves upward on the membrane by capillary action. In the case of a positive result the specific antibodies present on the membrane will react with the mixture conjugate and generate a coloured line. A green coloured line always appears in the control line and serves as verification that sufficient volume was added that proper flow was obtained and as an internal control for the reagents.

### PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed pouch until use.
- Do not use the test if pouch is damaged.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, do not eat, drink or smoke in the area.
- Clean up spills thoroughly using an appropriate disinfectant.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. Please, take note that the COVID-19 infection is very dangerous, extremely the precautions.
- The test should be discarded in a proper biohazard container after testing.
- The test must be carried out within 2 hours of opening the sealed bag.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control and test line zone) that are visible before using the test are completely normal. That does not mean failure on test functionality.
- Positive results should be processed following local laws and regulations.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch either at refrigerated or room temperature (2-30°C/36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

#### MATERIALS PROVIDED

- 20 Tests
- Instruction for use
- 1 Diluent (sample diluent)
- 20 Sterile Swabs
- 20 Plastic pipettes
- 20 Testing tubes or vials
- 1 SARS-CoV-2 control +
- 1 Laboratory rack

#### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

### SPECIMEN COLLECTION AND PREPARATION

Nasopharyngeal swab method:

- Bend shaft to follow curve of nasopharynx.
- Insert swab through nostril to posterior nasopharynx.
- Rotate swab a few times to obtain infected cells.
- For an optimal sample, repeat procedure using another nostril.

Send specimen to lab immediately (testing sensitivity decrease over time).

Cool specimen to 2°-8°C (36°-46.4°F) during storage and transport for 8 hours prior to testing.

### PROCEDURES


**Allow the tests, samples and buffer to reach room temperature (15-30°C/59-86°F) prior to testing. Do not open pouches until ready to perform the assay.**

**To process the collected nasopharyngeal swabs (see illustration 1):**

Use a separate testing tube or vial for each sample (swab). Add the diluent (15 drops) into the testing tube or vial, put the nasopharyngeal swab, mix and extract as much liquid possible from the swab. Remove the SARS-CoV-2 MonlabTest® from its sealed pouch and use it as soon as possible. Use a separate device for each sample. Dispense exactly 3 drops into the specimen well (S). Start the timer. Read the result at **10 minutes** after dispensing the sample.


### Illustration 1 Nasopharyngeal swab

### INTERPRETATION OF RESULTS

#### Illustration 2

**POSITIVE:** Two lines appear across the central window, a **red** test line marked with the letter T and a **green** control line marked with the letter C.

**NEGATIVE:** Only one **green** line appears across the control line region marked with the letter C (control line).

**INVALID:** Total absence of the **green** control coloured line regardless the appearance or not of the **red** test line. Note: Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor. See illustration 2.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured test line in the result line region (T) will vary depending on the concentration of antigens in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural controls are included in the test:

- A green line appearing in the control line region (C). It confirms sufficient specimen volume and correct procedural technique.

#### External Quality Control


Each kit contains a positive control material. Use the control swab to check that the extraction reagents and the test are working properly. Also use the controls to test that you are able to correctly perform the test procedure.

#### Quality Control Procedure:

**SARS-CoV-2 Positive control:** Remove the SARS-CoV-2 positive control from its sealed pouch. Add the diluent (15 drops) in a testing tube. Put the SARS-CoV-2 positive control swab, mix 60 seconds and extract as much liquid possible from the swab. Discard the swab. Remove the test from its sealed pouch and dispense 3 drops of the positive control liquid into the specimen well (S).

Result: SARS-CoV-2 positive (see interpretation of results).

### Illustration 3. Positive Control

### LIMITATIONS

1. SARS-CoV-2 MonlabTest® will only indicate the presence of Coronavirus (SARS-CoV-2) in the specimen (qualitative detection) and should be used for the detection of Coronavirus (SARS-CoV-2) antigens in nasopharyngeal specimens only (from swab). Neither the quantitative value nor the rate of increase in Coronavirus (SARS-CoV-2) antigens concentration can be determined by this test.
  2. Positive results do not rule out co-infections with other pathogens.
  3. VTM, UTM and Saline Buffer are the transport media validated for be used with this test (SARS-CoV-2 MonlabTest®). Always, following proportion 1:1 (transport media and sample diluent provided with the test). When using transport media the sensitivity of the test can be reduce due to excessive dilution of sample. The preference is to use the sample immediately after taking it.
1. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of Coronavirus (SARS-CoV-2) infection.
  2. This test provides a presumptive diagnosis of Coronavirus (SARS-CoV-2) respiratory infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

### EXPECTED VALUES

On March 11, 2020 the World Health Organization (WHO) declared COVID-19 a pandemic.

COVID-19 has led to an emerging, rapidly evolving situation: reaching more than 200 countries, areas, or territories, to date, and causing thousands of deaths, with an overall mortality around 3%.

Emerging data from patients afflicted by COVID-19 indicate that genetic immunologic, metabolic and environmental factors are involved in the pathogenesis of COVID-19. It was found that certain infected people were asymptomatic. These individuals are an important reservoir for spread of the infection.

### PERFORMANCE CHARACTERISTICS

#### Detection limit

The detection limit value (typical value) is: 1.0 ng/mL of SARS-CoV-2 Recombinant Nucleoprotein or  $1 \times 10^3$  TCID<sub>50</sub>/mL of 2019nCoV/USA-WA1/2020.

#### Sensitivity and Specificity

It was performed an evaluation with 262 nasopharyngeal samples from people suspected of infection by SARS-CoV-2.

The evaluation was performed using SARS-CoV-2 MonlabTest® versus qPCR technique.

IC test: SARS-CoV-2 MonlabTest®	qPCR technique			Total
	+	-		
+	26	1		27
-	2	233		235
Total	28	234		262

	SARS-CoV-2 MonlabTest® vs qPCR technique	
		95% CI (Confidence interval)
<b>Sensitivity</b>	92.9%	76.5-99.1%
<b>Specificity</b>	99.6%	97.6-100.0%
<b>PPV</b>	96.3%	81.0-99.9%
<b>NPV</b>	99.1%	97.0-99.9%

Recommendations of the WHO: "Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)" were considered. The sensitivity of the test was calculated with nasopharyngeal samples with high viral load in the range of Ag-RDT detection. High viral loads are expected in early symptomatic phases of the illness (first 5-7 days of illness).

### Hook Effect

SARS-CoV-2 MonlabTest® does not show hook effect at the concentration of SARS-CoV-2 protein tested. (202500.0 ng/mL)

### Cross-Reactivity

It was performed an evaluation to determine the cross reactivity of SARS-CoV-2 MonlabTest®. There is not cross reactivity with common respiratory pathogens, other organisms and other substances that could cause infections:

Adenovirus	Coronavirus strain OC43	Lactoferrin (human)	Salmonella enteritidis/paratyphi A/typhi/typhimurium
Astrovirus	Entamoeba histolytica	Legionella	Shigella boydii/dysenteriae/flexneri/sonnei
Calprotectin (human)	Escherichia coli O157:H7	Listeria monocytogenes	Streptococcus pneumococcal
Campylobacter jejuni	Giardia	Norovirus GI/Norovirus GII	Streptococcus pyogenes
Clostridium difficile GDH/Toxin A/Toxin B	Helicobacter pylori	Rotavirus	Transferrin (human)
Coronavirus strain 229E	Hemoglobin (human/bovine/pig)	Respiratory syncytial virus	Yersinia enterocolitica O:3/O:9
Coronavirus strain NL63	Influenza A/Influenza B		

SARS-CoV-2 MonlabTest® has some cross reaction with SARS and very low cross reactions with MERS.

### Interferences

It was performed an evaluation to determine the possible interferences of SARS-CoV-2 MonlabTest®. There are not interferences with the substances tested:

Metronidazole	Loratadine	Loperamide hydrochloride (Fortasec)	Phenoxymethylpenicillin potassium
Ampicillin	Dexchlorpheniramine (Polaramine)	Heparin (Hibor)	Ambroxol hydrochloride (Mucosan)
Osetamivir	Ebastine (Ebastel)	Almagato (Almax)	Macrogol 3350 (Movicol)
Amantadine	Acetyl Salicylic (Adiro)	Fosfomicin (Monurol)	Lysine Carbocysteinate (Pectox)
Ribavirin	Ibuprofen (Espidifen)	Acetylcysteine (Flumucil)	Hydroxyzine dihydrochloride
Codeine (Toseina)	Paracetamol (Dolocatil)	Dexketoprofen trometamol (Enantyum)	Lorazepam
Benzocaine (Angileptol)	Metamizole (Nolotil)	Levofloxacin	Amoxicillin
Cloperastine (Flutox)	Prednisone	Ciprofloxacin	Mercaptopurine
Carbocysteine (Miston mucolitica)	Omeprazole	Rifampicin (Rifaldin)	

### RECOMENDATIONS

Recommendations **World Health Organization**: Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 September 2020.

- To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.
- Where possible, all positive samples giving positive results should be transported to laboratories with NAAT capability for confirmatory testing.
- This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT-confirmed COVID-19 outbreaks.
- To monitor trends in disease incidence in communities.
- For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
- A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
- Even SARS-CoV-2 MonlabTest® test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so SARS-CoV-2 MonlabTest® could detect as positive.

### REFERENCES

- Palestino, G.; Garcia-Silva, I.; González-Ortega, O. and Rosales-Mendoza, S. **Can nanotechnology help in the fight against COVID-19? Expert Review of Anti-infective Therapy.** <https://doi.org/10.1080/14787210.2020.1776115>. Jun. 2020.
- Abiodun Iwalokun, B.; Olalekan, A.; Adenipekun, E.; Ojo, O.; Senapon Iwalokun, O.; Mutiu, B.; Orija, O.; Adegbola, R.; Salako, B. and Akinloye, O. **Improving the Understanding of Immunopathogenesis of Lymphopenia as a Correlate of SARS-COV-2 Infection Risk and Disease Progression in African Patients: UGLY SARS-COV-2 Study Protocol.** Journal of Medical Internet Research. <https://preprints.jmir.org/preprint/21242>. Jun. 2020.
- Giménez, E.; Albert, E.; Torres, I.; Remigia, M. J.; Alcaraz, M. J.; Galindo, M. J.; Blasco, M. L.; Solano, C.; Forner, M. J.; Redón, J.; Signes-Costa, J. and Navarro, D. **SARS-CoV-2-reactive interferon-γ-producing CD8+T cells in patients hospitalized with Coronavirus Disease 2019.** Journal of Medical Virology. <https://doi.org/10.1002/jmv.26213>. Jun. 2020.
- World Health Organization, Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays. Interim guidance. 11 September.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	Manufacturer		For <i>in vitro</i> diagnostic use only
	Authorized representative		Consult instructions for use
	Contains sufficient for <n> tests		Keep dry
	Catalogue Code		Temperature limitation
	Lot Number		Use by
	Diluent (sample diluent)		

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